


**medax**

**HEMAX**

## BONE MARROW BIOPSY SYSTEM WITH EXTRACTION CANNULA



Bone marrow biopsy system provided with special “HM” or “HF” safe and atraumatic extraction cannula.

- Coloured-code ergonomic handle.
- Available in different tip configuration: Diamond, Aspiration and Fish Mouth.
- “HM” or “HF” additional extraction cannula to get the sample in a completely atraumatic and safe way.
- Luer Lock connection for syringe.
- Remover guide provided with indicators to facilitate the sample expulsion and to check its length during the procedure. - Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

## BONE MARROW BIOPSY SYSTEM WITH EXTRACTION CANNULA

GAUGE	COLOUR	LENGTH (MM)		
8	● Dark Green	08 070	08 100	08 150
9	● Pink	09 070	09 100	09 150
11	○ White	11 070	11 100	11 150
13	● Light blue	13 070	13 100	13 150

Different sizes on demand

### Ordering information:

HEMAX "M" DIAMOND TIP	H	M	GAUGE	LENGTH	0	0
HEMAX "M" BONE MARROW ASPIRATION	H	M	GAUGE	LENGTH	0	1
HEMAX "M" FISH MOUTH TIP	H	M	GAUGE	LENGTH	0	2
HEMAX "F" DIAMOND TIP	H	F	GAUGE	LENGTH	0	0
HEMAX "F" BONE MARROW ASPIRATION	H	F	GAUGE	LENGTH	0	1
HEMAX "F" FISH MOUTH TIP	H	F	GAUGE	LENGTH	0	2

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

