

medax

**MED-L**

BONE MARROW BIOPSY AND ASPIRATION SYSTEM  
WITH SPECIAL EXTRACTION CANNULA



EXTRACTION CANNULA  
"P" VERSION

EXTRACTION CANNULA  
"S" VERSION

Bone marrow biopsy and aspiration system provided with special safe and atraumatic extraction cannula.

- Coloured-code ergonomic handle.
- Additional extraction cannula to get the sample in a completely atraumatic and safe way.
- "S" version optimal for multiple procedures.
- Luer Lock connection for syringe.
- Remover guide provided with indicators to facilitate the sample expulsion and to check its length during the procedure.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

# MED-L

## BONE MARROW BIOPSY AND ASPIRATION SYSTEM WITH SPECIAL EXTRACTION CANNULA

GAUGE	COLOUR	LENGTH (MM)		
8	● Dark Green	08 100	08 150	
11	○ White	11 100	11 150	
13	● Light blue	13 070	13 100	13 150

Med-L "P" version 13G not available

### Ordering information:

MED-L "S" VERSION DIAMOND TIP

X	L	GAUGE		LENGTH		0	S
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

MED-L "P" VERSION FISH MOUTH TIP

X	L	GAUGE		LENGTH		0	P
		<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>		

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744 • email: [customercare@medax.it](mailto:customercare@medax.it) • PEC: [medax@legalmail.it](mailto:medax@legalmail.it) • [www.medax.it](http://www.medax.it)

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.