

medax

HEPAMED

MENGHINI KIT FOR LIVER BIOPSY



Histological aspiration biopsy kit using the Menghini technique without ultrasound imaging.

- Syringe of 12ml with self-locking system to guarantee a constant vacuum.
- Ultra thin sharpened cannula with high performance steel for a controlled penetration and histological aspiration.
- Scalpel to facilitate the percutaneous needle insertion.
- Hypodermic needle.
- Centimeter marked cannula and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

HEPAMED

MENGHINI KIT FOR LIVER BIOPSY

| GAUGE | COLOUR | LENGTH (MM) |
|-------|---------|-------------|
| 15 | ● Blue | 15 088 |
| 16 | ○ White | 16 088 |
| 17 | ● RED | 17 088 |
| 18 | ● Pink | 18 088 |

Ordering information:

HEPAMED



Single box: 10 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744 • email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.