

medax

SONOMED

FULL-CORE BIOPSY SYSTEM
FOR HISTOLOGICAL ASPIRATION BIOPSY



SONOMED, full-core biopsy system for ultrasound guided histological aspiration biopsy.

- Ultra thin sharpened cannula with trocar stylet for a quick and painless percutaneous insertion.
- Syringe of 12ml.
- Compatible with SONOGUN.
- Sterilized by ETO, shelf life 5 years.

medax

SONOGUN

REUSABLE ASPIRATION GUN



SONOGUN, Reusable aspiration gun usable with disposable needle SONOMED.

- Light Weight
- One hand usage.
- Compatible with 10ml and 12ml syringe.





medax
medical devices

Your partner in Biopsy and Special Needles

medax

SONOMED

FULL-CORE BIOPSY SYSTEM FOR HISTOLOGICAL ASPIRATION BIOPSY

GAUGE	COLOUR	LENGTH (MM)	
17	 RED	17 100	17 160
18	 Pink	18 100	18 160
20	 Yellow	20 100	20 160
21	 Light green	21 100	21 160

Ordering information:

SONOMED



Single box: 10 PCS.

SONOGUN



Single box: 1 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744 • email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.