

medax

VILLOX

**SPECIAL SHARPENED NEEDLES FOR TRANSABDOMINAL
RETRIEVAL OF CHORIONIC VILLI DURING PRENATAL
DIAGNOSIS**



Special sharpened needles for transabdominal retrieval of chorionic villi during prenatal diagnosis.

- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Internal or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.


medax
medical devices

Your partner in Biopsy and Special Needles

medax

VILLOX

SPECIAL SHARPENED NEEDLES FOR TRANSABDOMINAL RETRIEVAL OF CHORIONIC VILLI DURING PRENATAL DIAGNOSIS

GAUGE	COLOUR	LENGTH (MM)		
19	 Light Brown	19 160	19 180	19 200
20	 Yellow	20 160	20 180	20 200

Ordering information:

VILLOX



Single box: 20 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.
In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.