


**medax**

# MEDEASY

SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM








Semi-automatic spring loaded biopsy system with selectable penetration depth.

- Adjustable penetration depth (10/20 mm).
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Sword tip.
- Can be used with CE introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

## SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)									
12	 Light blue	12 060	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300	
14	 Green	14 060	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300	
16	 White	16 060	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300	
18	 Pink	18 060	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300	
20	 Yellow	20 060	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300	

Different sizes on demand - 60mm length not available in configuration of coaxial needle and blunt tip stylet.

### Ordering information:

MEDEASY BEVEL TIP	M S	GAUGE	LENGTH	0 0
MEDEASY BEVEL TIP "LIGHT VERSION"	M L	GAUGE	LENGTH	0 0
MEDEASY SWORD TIP	M S	GAUGE	LENGTH	0 2
MEDEASY SWORD TIP "LIGHT VERSION"	M L	GAUGE	LENGTH	0 2
INTRODUCER COAXIAL NEEDLE FOR MEDEASY	C E	GAUGE	LENGTH	0 0
MEDEASY BEVEL TIP WITH COAXIAL NEEDLE	M S	GAUGE	LENGTH	C 0
MEDEASY BEVEL TIP "LIGHT VERSION" WITH COAXIAL NEEDLE	M L	GAUGE	LENGTH	C 0
MEDEASY BEVEL TIP WITH COAXIAL NEEDLE INCLUDING BLUNT TIP STYLET	M S	GAUGE	LENGTH	K 0
MEDEASY BEVEL TIP "LIGHT VERSION" WITH COAXIAL NEEDLE INCLUDING BLUNT TIP STYLET	M L	GAUGE	LENGTH	K 0

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

