

medax

UNIMED

SOFT TISSUE BIOPSY NEEDLE

Soft tissue biopsy needles compatible with different reusable biopsy guns.

- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CU introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

GAUGE	COLOUR	LENGTH (MM)							
12	● Light blu	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand - UN, UP & UF not available in 12 Gauge.



COMPATIBLE WITH

- M-GUN
- PROMAG 2.2™
- ABS

Biopsy needle with Beveled tip

U	N	GAUGE	LENGTH	C	1
				0	0

Biopsy needle with Trocar tip

U	N	GAUGE	LENGTH	C	1
				0	1

Biopsy needle with Beveled tip including Introducer coaxial

U	N	GAUGE	LENGTH	C	1
				C	0

Biopsy needle with Trocar tip including Introducer coaxial

U	N	GAUGE	LENGTH	C	1
				C	1



COMPATIBLE WITH

- MEDGUN
- MAGNUM®

Biopsy needle with Beveled tip

U	P	GAUGE	LENGTH	C	1
				0	0

Biopsy needle with Trocar tip

U	P	GAUGE	LENGTH	C	1
				0	1

Biopsy needle with Beveled tip including Introducer coaxial

U	P	GAUGE	LENGTH	C	1
				C	0

Biopsy needle with Trocar tip including Introducer coaxial

U	P	GAUGE	LENGTH	C	1
				C	1

medax
medical devices

Your partner in Biopsy and Special Needles

UNIMED

SOFT TISSUE BIOPSY NEEDLE



COMPATIBLE WITH
 • MEDGUN
 • MAGNUM®

Biopsy needle with Beveled tip	U	D	GAUGE	LENGTH	0	0
Biopsy needle with Trocar tip	U	D	GAUGE	LENGTH	0	1
Biopsy needle with Beveled tip including Introducer coaxial	U	D	GAUGE	LENGTH	C	0
Biopsy needle with Trocar tip including Introducer coaxial	U	D	GAUGE	LENGTH	C	1



COMPATIBLE WITH
 • MEDGUN
 • MAGNUM®

Biopsy needle with Beveled tip	U	G	GAUGE	LENGTH	0	0
Biopsy needle with Trocar tip	U	G	GAUGE	LENGTH	0	1
Biopsy needle with Beveled tip including Introducer coaxial	U	G	GAUGE	LENGTH	C	0
Biopsy needle with Trocar tip including Introducer coaxial	U	G	GAUGE	LENGTH	C	1



COMPATIBLE WITH
 • MEDGUN
 • MAGNUM®

Biopsy needle with Beveled tip	U	G	GAUGE	LENGTH	S	0
Biopsy needle with Trocar tip	U	G	GAUGE	LENGTH	S	1
Biopsy needle with Beveled tip including Introducer coaxial	U	G	GAUGE	LENGTH	E	0
Biopsy needle with Trocar tip including Introducer coaxial	U	G	GAUGE	LENGTH	E	1



COMPATIBLE WITH
 • MEDCORE
 • MAGNUM®

Biopsy needle with Beveled tip	U	C	GAUGE	LENGTH	0	0
Biopsy needle with Trocar tip	U	C	GAUGE	LENGTH	0	1
Biopsy needle with Beveled tip including Introducer coaxial	U	C	GAUGE	LENGTH	C	0
Biopsy needle with Trocar tip including Introducer coaxial	U	C	GAUGE	LENGTH	C	1

UNIMED

SOFT TISSUE BIOPSY NEEDLE



COMPATIBLE WITH

- MEDCORE
- MEDGUN
- MAGNUM®
- PALLIUM

Biopsy needle with Beveled tip

U	C	GAUGE	LENGTH	S	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip

U	C	GAUGE	LENGTH	S	1
---	---	-------	--------	---	---

Biopsy needle with Beveled tip including Introducer coaxial

U	C	GAUGE	LENGTH	E	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip including Introducer coaxial

U	C	GAUGE	LENGTH	E	1
---	---	-------	--------	---	---



COMPATIBLE WITH

- MEDCORE
- MEDGUN

Biopsy needle with Beveled tip

U	C	GAUGE	LENGTH	G	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip

U	C	GAUGE	LENGTH	G	1
---	---	-------	--------	---	---

Biopsy needle with Beveled tip including Introducer coaxial

U	C	GAUGE	LENGTH	F	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip including Introducer coaxial

U	C	GAUGE	LENGTH	F	1
---	---	-------	--------	---	---



COMPATIBLE WITH

- MEDMAG
- PROMAG™ 2.5 ULTRA

Biopsy needle with Beveled tip

U	M	GAUGE	LENGTH	0	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip

U	M	GAUGE	LENGTH	0	1
---	---	-------	--------	---	---

Biopsy needle with Beveled tip including Introducer coaxial

U	M	GAUGE	LENGTH	C	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip including Introducer coaxial

U	M	GAUGE	LENGTH	C	1
---	---	-------	--------	---	---



COMPATIBLE WITH

- EVO CORE®
- CORAZOR®
- MAGNUM®

Biopsy needle with Beveled tip

U	T	GAUGE	LENGTH	0	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip

U	T	GAUGE	LENGTH	0	1
---	---	-------	--------	---	---

Biopsy needle with Beveled tip including Introducer coaxial

U	T	GAUGE	LENGTH	C	0
---	---	-------	--------	---	---

Biopsy needle with Trocar tip including Introducer coaxial

U	T	GAUGE	LENGTH	C	1
---	---	-------	--------	---	---

UNIMED

SOFT TISSUE BIOPSY NEEDLE



COMPATIBLE WITH

- PALLIUM
- MAGNUM®
- MEDCORE
- MEDGUN

Biopsy needle with Beveled tip

U	L	GAUGE		LENGTH		0	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip

U	L	GAUGE		LENGTH		0	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Beveled tip including Introducer coaxial

U	L	GAUGE		LENGTH		C	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip including Introducer coaxial

U	L	GAUGE		LENGTH		C	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		



COMPATIBLE WITH
• FAST GUN® LIGHT

Biopsy needle with Beveled tip

U	F	GAUGE		LENGTH		0	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip

U	F	GAUGE		LENGTH		0	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Beveled tip including Introducer coaxial

U	F	GAUGE		LENGTH		C	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip including Introducer coaxial

U	F	GAUGE		LENGTH		C	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		



COMPATIBLE WITH
• DANA 2.2MG

Biopsy needle with Beveled tip

U	H	GAUGE		LENGTH		0	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip

U	H	GAUGE		LENGTH		0	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Beveled tip including Introducer coaxial

U	H	GAUGE		LENGTH		C	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip including Introducer coaxial

U	H	GAUGE		LENGTH		C	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		



COMPATIBLE WITH
• MEDMAG
• PROMAG™ 2.5 ULTRA

Biopsy needle with Beveled tip

U	O	GAUGE		LENGTH		0	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip

U	O	GAUGE		LENGTH		0	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Beveled tip including Introducer coaxial

U	O	GAUGE		LENGTH		C	0
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Biopsy needle with Trocar tip including Introducer coaxial

U	O	GAUGE		LENGTH		C	1
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744
email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

