

**medax**

# UNIVERSAL

M-GUN COMPATIBLE BIOPSY NEEDLE  
ALSO USABLE WITH MOST COMMON MARKET  
REUSABLE BIOPSY SYSTEMS



M-Gun compatible needle also usable with the most common reusable guns:

- Angiotech/Pro-mag™ 2.2
- Bard®/Biopty®
- Gallini/ABS
- US Biopsy/MBD-23
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Trocar tip.
- Can be used with CU introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.

**medax**  
medical devices

Your partner in Biopsy and Special Needles

# UNIVERSAL

M-GUN COMPATIBLE BIOPSY NEEDLE ALSO USABLE WITH AND MOST COMMON MARKET REUSABLE BIOPSY SYSTEMS

GAUGE	COLOUR	LENGTH (MM)								
		080	100	120	160	180	200	250	300	
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300	
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300	
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300	
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300	

Different sizes on demand

### Ordering information:

UNIVERSAL NORMAL TIP

U	N	GAUGE		LENGTH			0	0

UNIVERSAL TROCAR TIP

U	N	GAUGE		LENGTH			0	1

INTRODUCER COAXIAL NEEDLE FOR UNIVERSAL

C	U	GAUGE		LENGTH			0	0

UNIVERSAL WITH COAXIAL NEEDLE IN THE SAME PACKAGE

U	N	GAUGE		LENGTH			C	0

Single box: 10 PCS.

**Medax Srl Unipersonale**

**Headquarters:** Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy  
 Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744  
 email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

**Registered Office:** Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy  
 Vat N./Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.

