


medax

MEDWIRE

BREAST LOCALIZATION NEEDLE
FOR NON PALPABLE MAMMARY LESIONS



SINGLE HOOK

MARKED VERSION

STIFFENED VERSION

Breast localization needle for non palpable mammary lesions.

- Available in SINGLE wire configuration.
- Also available in marked and stiffened versions.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.

medax
medical devices

Your partner in Biopsy and Special Needles

medax

MEDWIRE

BREAST LOCALIZATION NEEDLE FOR NON PALPABLE MAMMARY LESIONS

VERSION HOOK	GAUGE Ø	WIRE (MM)	LENGTH (MM)					
L-Single	20	0,25-0,30	20 060	20 080	20 100	20 110	20 120	20 160
L-Single	21	0,25	21 060	21 080	21 100	21 110	21 120	21 160

Ordering information:

MEDWIRE SINGLE HOOK WIRE

X	W	GAUGE	LENGTH	0	S
---	---	-------	--------	---	---

MEDWIRE SINGLE HOOK WIRE MARKED VERSION

X	W	GAUGE	LENGTH	1	S
---	---	-------	--------	---	---

MEDWIRE SINGLE HOOK WIRE STIFFENED VERSION

X	W	GAUGE	LENGTH	2	S
---	---	-------	--------	---	---

Single box: 20 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744 • email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates. All products undergo intensive clinical testing and are fully EC and FDA approved.