

medax

AMNIOX

PRENATAL DIAGNOSIS NEEDLE



Special sharpened needles for transabdominal retrieval of amniotic fluid and chorionic villi during prenatal diagnosis.

- Transparent Luer Lock connector.
- Sharpened cannula for a safe and painless insertion and penetration.
- Internal or external echogenic treatment for a correct positioning under ultrasound guidance.
- Centimeter markings and sliding stopper.
- Sterilized by ETO, shelf life 5 years.







medax
medical devices

Your partner in Biopsy and Special Needles

medax

AMNIOX

PRENATAL DIAGNOSIS NEEDLE

GAUGE	COLOUR	LENGTH (MM)		
18	 Pink	18 100	18 150	18 200
19	 Light brown	19 100	19 150	19 200
20	 Yellow	20 100	20 150	20 200
21	 Light green	21 100	21 150	21 200
22	 Black	22 100	22 150	22 200
23	 Light Blue	23 100	23 150	23 200

Ordering information:

AMNIOX

A	M	GAUGE		LENGTH			0	0

Single box: 20 PCS.

Medax Srl Unipersonale

Headquarters: Via S. Pertini, 4 • 41039 • San Possidonio (MO) • Italy
Company direct No. : +39 0535 1812757 • Fax No : +39 0535 1812744 • email: customercare@medax.it • PEC: medax@legalmail.it • www.medax.it

Registered Office: Via R. Piva, 1/A • 46025 • Poggio Rusco (MN) • Italy
Vat N. /Fiscal Code N. Iscriz. Reg. Impr.: MN 02669860369 • N. REA: MN 233527, MO 403036 • Capitale Sociale Euro 100.011,00 i.v.



Medax's quality management system is certified to ISO 13485: 2016 standards.



In accordance with the requirements of the medical device 93/42/EEC directive and its relevant updates.
All products undergo intensive clinical testing and are fully EC and FDA approved.